SUPPORTING ROLL WITH MASSAGE

English, Français

CONGRATULATIONS

on your purchase of an ObusForme Supporting Roll with Massage. The versatile ObusForme Supporting Roll with Massage provides excellent support for the curvature of your lower back or neck, and can provide relief from pain aggravated by poor posture including back pain, neck pain, shoulder tension and headaches. Small, lightweight and portable, the ObusForme Supporting Roll with Massage is perfect for use anywhere you sit – at home, in the office and on the go!

FEATURES OF YOUR OBUSFORME SUPPORTING ROLL WITH MASSAGE

- Provides excellent, versatile support for the curvature of your lower back or neck
- Provides a choice of low or high intensity massage for the curvature of your lower back or neck
- Can provide relief from pain aggravated by poor posture including back pain, neck pain, shoulder tension and headaches

USING YOUR OBUSFORME SUPPORTING ROLL WITH MASSAGE

- To activate the massage option, adjust the switch to (I) low intensity or (II) high intensity depending on your preferred setting
- Your ObusForme Supporting Roll with Massage can support the natural curvature of your neck while you sleep on your back
- A convenient elastic strap holds your ObusForme Supporting Roll with Massage in the position you desire; simply adjust and secure it using the button clasps

CARING FOR YOUR OBUSFORME SUPPORTING ROLL WITH MASSAGE

- Do not allow the interior foam or the battery compartment of your Supporting Roll with Massage to become damp or wet.
- Do not rub the cover excessively

BATTERY REPLACEMENT

1. Locate the battery compartment on the back of the Supporting Roll with Massage.
2. Remove the exhausted batteries from the Supporting Roll with Massage.
3. Insert four, new AA (LR6) alkaline batteries as indicated inside the battery compartment.

Helpful Hints:

- For longer battery life, we recommend use of alkaline batteries.
- Fit the tab on the battery compartment cover into the slot in the foam and close the cover.
- Do not mix old and new batteries.
- Do not mix different types of batteries: alkaline, standard (carbon-zinc) or rechargeable (nickel-cadmium).
- Remove the batteries during long periods of non-use. Always remove exhausted batteries from the Supporting Roll with Massage. Battery leakage and corrosion can cause damage. Dispose of batteries safely.
- Never short circuit the battery terminals.
- Non-rechargeable batteries are not to be recharged.
- Only batteries of the same or equivalent type as recommended in the "Battery Replacement" instructions are to be used.
- If removable rechargeable batteries are used, they are only to be charged under adult supervision.
- Rechargeable batteries are to be removed from the Supporting Roll with Massage before they are charged.
- If you use a battery charger, it should be examined regularly for damage to the cord, plug, enclosure and other parts. Do not use a damaged battery charger until it has been properly repaired.

OBUSFORME LIMITED WARRANTY

ObusForme guarantees all items are free from defects in workmanship and materials for a period of time between the original purchase date and that stated below. This guarantee applies when items are used for the purpose intended. Items will be repaired/replaced (at our option), with new, refurbished parts/products and/or substitutes, if the ORIGINAL purchaser has completed and returned the Warranty Registration within 30 days of purchase and includes original receipt. Shipping, customs, duties and taxes must be PRE-PAID TO and FROM ObusForme by the PURCHASER. This warranty gives you rights that vary by province/state. This warranty may change.

WHAT IS NOT COVERED

Wear and tear, aging (including foam discoloration, flattening, density, consistency), accidental damages, alterations, mishandling, faulty adjustment, misuse, improper care, power damage, rental use, discontinued items, service by anyone other than ObusForme, shipping damages, over inflation, neglect, items sold 'as is' or damage due to natural acts are NOT covered.

WARRANTY TIME FRAME

Backrest Frame: Lifetime

(Cover, foam, lumbar pad and other Backrest parts/materials are NOT covered)

Seat Frame: 1 year

(Cover, foam and other Seat parts/materials are NOT covered)

Back Therapy: 1 year

(Back Therapy includes Backlife, Back Belts, Back Packs, Drivers Seats)

Sleep/Foot/Muscle Therapy: 1 year

(Pillow cases are NOT covered)

Electrical Parts: 1 year

(This includes wires, adaptors, plugs and other electrical parts/components)

HOW TO OBTAIN WARRANTY SERVICE

Please obtain a Return Authorization Number and instructions prior to sending your item or it may be denied. Please inform Customer Service by:

Mail:

HoMedics Group Canada

344 Consumers Road

Toronto, Ontario, Canada M2J 1P8

Tel: (416) 785-1386 Fax: (416) 785-5862

Toll Free: 1-888-225-7378

8:30 a.m. to 5:00 p.m., Mon - Fri EST

www.obusforme.com

OBUSFORME®
THE BODYCARE EXPERTS™

ASSOCIATION OF CANADIAN ERGONOMISTS

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