Uplift Commode Assist:
The lifting commode chair that gives independence and privacy when it's needed most.

Uplift Commode Assist is a selfpowered, lifting commode chair that helps users gently seat and raise themselves. It provides safe, controlled support and lifts up to 70% of a user's weight.

Versatile and sturdy, the Commode Assist can be used standalone or over the toilet. For safety, it automatically locks into position when fully depressed.

Benefits

- Helps keep people active and independent.
- Relieves stress and strain on joints.
- · Lifts only as needed, allowing users to maintain muscular strength.

Who can benefit

- People who have lost sufficient upper or lower body strength to get up unassisted.
- Those who have lost mobility as a result of:
- Muscular dystrophy
- Lou Gehrig's disease (ALS)
- Post-operative conditions
- Degenerative joint disease
- Arthritis
- Parkinson's disease









The Commode Assist eases you gently into a seated position and automatically locks in place. To stand, simply press down on the seat release lever (see below).

Using both hands to support yourself, simply shift forward and begin to stand.

The lifting mechanism activates automatically and provides a gentle and stable lift up from the commode.

Commode Assist Specifications Base Dimensions Depth:

Seat Height

21"/53 cm

Width: 26"/66 cm 16.75"/43 cm

Seat Dimensions Depth: Width: 16.5"/42 cm

Minimum: 17.75"/45 cm

Maximum: 23.5"/60 cm

Lift 9"/23 cm @ 20° Weight 21 lb/9.5 kg

Easy to clean snap-on seat. Self-powered lifting mechanism releases slowly as user stands. Extendable legs with five height settings. Non-skid rubber tips for added safety. Steel frame is white powder coated to protect against moisture. Model No. Description Weight Range

Bucket slides out for easy cleaning. Includes splash guard and cover.

Cushioning pads enhance comfort.

Removable back

easy positioning over toilet.

rest allows for

Seat release lever

commode assist

locks when

depressed.

& auto-lock hook -



simply push down on the seat release lever. Now both hands are free to use for support when rising (see above sequence).

CA 200* Lifting commode

80-300+ lb chair 35-136+ kg

*some assembly required



Bringing quality to life.