The Obus Forme® Soft Gel Hot & Cold Compress

This easy-to-use, reusable compress helps provide soothing relief for strained and aching muscles, head, neck, tooth and back aches, bruises, abrasions, minor burns, insect bites and other injuries. Please read this carefully before use.

Cold Pain Relief

Recommended directly after an injury and each day after until the swelling and pain subside. Using the cold compress will lessen damage to the tissue, decrease swelling and consequently, help reduce pain. Following an injury, gradually return to your former activity level.

For optimum effectiveness of cold therapy, this compress should be stored in the freezer for a minimum of 2 hours. This compress will not freeze solid; it remains pliable to conform to the shape of the body.

When applying the compress, wrap it in a thin cloth that will not impede cold transfer. Cold therapy can be applied for 5-20 minutes at a time, depending on the location and size of the injury.

Hot Pain Relief

Not recommended for at least 2 days after an injury. Heat allows the pain to subside temporarily, helping to warm and loosen your joints so you can feel more comfortable moving your injured body part.

Stovetop: Immerse the compress in hot or boiling water for 5 to 7 minutes.

Microwave: Fully immerse the compress in water in a microwave-safe container. Heat on high power for 3 minutes. If further heating is required, microwave at 30-second intervals.

Safely remove the heated compress from the water and wrap it in a thin cloth that will not impede heat transfer. Hot therapy can be applied for as long as 20-30 minutes. Comfort level and heat tolerance should be used as guidelines.

CAUTION!

If you have concerns regarding your health, please consult your physician before use. An injury that is really painful or pain that does not resolve itself after 24-48 hours of self-care should be checked by a professional. Do not apply the compress directly to the skin. To avoid burns, use the hot compress carefully. Use only as directed. Individuals with circulatory or diabetic problems should consult a physician before using this product. Do not ingest. Do not puncture. Contents may stain. Non-toxic – Safe for use with children.

For additional information on back pain and general pain relief, new products and special offers or to subscribe to our newsletter, please visit www.obusforme.com

Corporate Member of ACE
Membre corporatif de l'ACE

ASSOCIATION OF CANADIAN ERGONOMISTS
ASSOCIATION CANADIENNE D'ERGONOMIE

HoMedics Group Canada, 344 Consumers Road, Toronto, Ontario M2J 1P8

Tel: (416) 785-1386 Fax: (416) 785-5862 Toll Free: 1-888-225-7378 www.obusforme.com

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