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
Rev 1 1008

## Digital Scale


## User's manual

Model No. 116-970

Rev 1 1008

## Your Scale at a Glance

- 1 - LCD
- 2 - Front electrodes
- 3 - Back electrodes
- 4 - "Up" button
- 5 - "Set" button
- 6 - "Down" button

## Precautions


- Do not use this scale if you are pregnant or if you have a pacemaker or any other medical implant. The current passing through your body may cause your medical equipment to malfunction.
- This scale is for home use only and should not be used as a medical device.
- To reduce the risk of falls, do not stand at the edge of the scale. Do not use this scale on a floor that is slippery or wet.
- Persons who are disabled should be helped when using this scale.

**Note:** Body fat readings vary with body water percentages. The following factors will also affect body fat percentage readings: diet, alcohol consumption, heavy exercising, certain prescription drugs or diuretics, illness and menstrual cycle.

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## Precautions

- Before use, remove your socks and ensure your feet are clean and dry.
- Do not wear clothes when taking a reading, as clothing could influence reading accuracy.
- Over-dieting and/or water loss, as well as intense exercise, will influence readings negatively.
- Proper contact with the metal surfaces (electrodes) is necessary to obtain a proper reading. Keep the scale clean and only step on it with clean, bare feet.
- Do not move while the scale is taking measurements.
- This scale will not provide accurate results for people who are younger than 10 years, people who are very muscular or athletic (i.e. body builders) pregnant women and people with dropsy.

## Instructions for set-up

### Before using the scale

1. Set the switch in the back of the scale to measure in kg (kilograms) or lb (pounds).
2. The battery compartment is located on the back of the scale. Before first use, pull out the insulation tape to allow battery contact. When the LCD indicates "LO", open the battery compartment and replace the battery with a new one. **Attention:** the battery polarity needs to be properly installed; otherwise, the scale will not work.
3. Put the scale on a hard, flat surface for accurate operation and safety. To avoid falling, do not put your foot on the edge of the scale.

**Note:** Before using the scale, select one personal number (P0 – P9) and input your specific parameters (gender, height, age). Without this information, the scale will only be able to detect your weight, but not your correct body fat and water percentages.

You will only need to enter this information once; the scale will store it under your personal number. Remember it and use it everytime you use the scale.

## Customizing and Using Your Scale

1. Press the "up" key to enter the setting mode.

2. Press the "up" or "down" key to reach the desired number (0 – 9). Press "set" to confirm.


3. Press the "up" or "down" key to select gender. Press "set" to confirm.

4. Repeat the steps above to set height and age. The parameters will then flash and be recorded.

5. When "0.0" appears, step on the scale, with your feet evenly on the 4 electrodes, as illustrated.

6. The scale first displays your weight.

7. 0000 indicates that the scale is analyzing (stay on the scale).

8. The scale then displays, in alternance, your weight, body fat, body water (TBW) and muscle mass (MUS) 4 times before shutting off automatically.

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## Understanding Body Fat, Water and Muscle Percentages

Due to differences in race and region, the chart below is for reference purposes only. **Attention:** in order to obtain more accurate readings, stand on the scale with bare feet, ensuring that your feet are clean, dry, and in good contact with the electrodes. Flashing lines under your body fat percentage will indicate if you are underfat (-), healthy (0), overfat (+) or obese (++).

### What is body fat percentage?

It is simply the percentage of fat that your body contains. If you are 150 pounds heavy and 10% fat, it means that your body contains 15 pounds of fat and 135 pounds of lean body mass (bone, muscle, organs, blood and everything else).

Knowing your body fat percentage can help you determine if your weight loss goals are realistic. Remember, weight loss does not always mean fat loss.

### How does the Physio Logic® Digital Scale measure body fat?

While body fat can be measured in many ways, the method used by this scale is BIA (bioelectrical impedance analysis). Body impedance is measured when a small amount of electric current is sent from the footpads through the body. The electric current flows more quickly through

water and muscle than it will through bone and fat. The scale measures the speed of this current. Based on this number, and along with a person's height, weight, age and gender, it is possible to calculate the percentage of body fat and water.

### Factors that may affect body fat percentage

Your body fat percentage may be affected by body hydration and by how long you have been awake / active. To ensure an accurate reading, measure your body fat at night. As body fat readings vary

### Body fat % norms for women

Age	20 – 39	40 – 59	60 – 79
Low fat	< 21 %	< 23 %	< 24 %
Average	21 – 33 %	23 – 35 %	24 – 36 %
High fat	33 – 39 %	35 – 40 %	36 – 42 %
Obese	> 39 %	> 40 %	> 42 %

### Body fat % norms for men


Age	20 – 39	40 – 59	60 – 79
Low fat	< 8 %	< 11 %	< 13 %
Average	8 – 19 %	11 – 22 %	13 – 25 %
High fat	19 – 25 %	22 – 28 %	25 – 30 %
Obese	> 25 %	> 28 %	> 30 %

# Understanding Body Fat, Water and Muscle Percentages

throughout the day, measuring at roughly the same time of the day and under consistent conditions will yield the best results.

The following factors also affect body fat percentage: diet, alcohol consumption, intense exercise, certain prescription drugs or diuretics, illness, menstrual cycle.

## Example of body fat fluctuation throughout the day

## What is the Body Mass Index?

The body mass index (BMI) is a number calculated from a person's **weight** and **height**. The BMI provides a reliable indicator of body fat for most people and is used to screen for weight categories that may lead to health problems.

Body mass index calculation is very straightforward and only requires two measurements: your height and your weight.

## How to calculate your BMI

Metric units:

$$BMI = \frac{\text{Weight (kg)}}{\text{Height (m)} \times \text{Height (m)}}$$

Example 1: Someone who is 1.70 m tall and weighs 60 kg:  $60 / (1.7 \times 1.7) = 20.8$  their BMI of 20.8 is in the Normal category.

English units:

$$BMI = \frac{\text{Weight (lb)}}{\text{Height (in)} \times \text{Height (in)}} \times 703$$

Example 2: Someone who is 5'6" tall (5'6" = 66 inches) and weighs 160 lb:  $(160 / (66 \times 66)) \times 703 = 25.8$  this person is in the overweight category.

BMI	Weight Status
< 18.5 %	Underweight
18.5% - 24.9 %	Normal
25% - 29.9 %	Overweight
> 30 %	Obese

## What is Body Water?

Body water is all of the water content in the human body. Lean muscle tissue is approx. 75% water. Blood contains 95% water, body fat 14% water and bone 22%.

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
## Cleaning and maintenance

- Always place the scale on a clean and level surface.
- Do not use this scale in very high or low temperatures.
- Do not jump or stomp on the scale, and avoid dropping it as it is a highly precise electronic instrument.
- To prevent scratches on the glass panel, do not take measurements while wearing any type of footwear.
- Clean the scale using a soft cloth. Avoid the use of any harsh chemicals so as to not damage the scale.
- This scale is for home use only. Do not use in a commercial setting. The body fat % chart is supplied for reference purposes only.
- Remove the battery from the scale if you plan to leave it unused for a certain period of time.

## Troubleshooting

### The scale powers off after only a few seconds.

- Check the battery: ensure that it is new and installed correctly.
- Ensure that your feet are in proper contact with the metal electrodes.

### The screen displays “----”.

- The scale is setting its zero value; please wait until “0” is displayed.
- No personal parameters were found; please customize your scale (see page 4).

### Only my weight is displayed (no body fat percentage).

- Shaking / movement will disturb body fat measurement; keep your body straight when weighing, and ensure that your bare feet are in good contact with the metal plates (electrodes).

### The indicated body fat percentage is much higher than normal.

- Verify the personal parameters entered (gender, age and height). Ensure that your feet are bare and dry, and that they are in good contact with the electrodes.
- Body fat readings may be affected by a major water loss, a state of agitation or callused feet.

## 5 Year Limited Warranty

This product is warranted against defects in materials and workmanship for a period of 5 years. This warranty does not cover product abuse or normal wear and is for the original user only.

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