



Product Design Advantages



Support exactly where you need it
The 3D flat-knitted fabric moulds to natural body contours – creating specific compression zones.



No constriction
Elastic edges prevent constriction.



Technogel® Pads – Gentle massage for affected areas
Extremely soft, elastic material allows user to massage affected area.



Antibacterial – odour-inhibiting – skin-friendly
Anti-bacterial fibres integrated into fabric creating excellent wearing comfort. Inhibits the growth of bacteria and noticeably reduces odour formation.



Breathable
Elasticity of fibres absorb moisture quickly and transport it to the surface of the flat-knitted fabric, where it evaporates quickly.

Our supports are free of natural latex and can be machine washed at 30°. The antibacterial effect is retained even after multiple washing cycles.



Sensa Line Supports

Right for any situation:
Protecting the body while staying mobile, requires targeted support. Innovative Sensa Line supports help you stay active and mobile. They let you move confidently and securely under strain without restrictions; and gently support the healing process by helping to reduce pain.

Right for all parts of the body:
Whether you require support for your hand, knee, ankle, spine or elbow, Sensa Line has the right support for you. Sensa Line fits perfectly and is very skin-friendly. Antibacterial fibres integrated into the knitted fabric effectively inhibit the growth of bacteria to reduce odour formation.

With the Sensa Line supports, you get the best possible fit for every situation and every part of the body.

Sensa Line. Fits me perfectly.



Sensa Line



Ottobock®
QUALITY FOR LIFE

5470 Harvester Road, Burlington, ON L7L 5N5
Phone 800.665.3327 • Fax 800.463.3659
www.ottobock.ca • info@ottobockcanada.com

*All trademarks in this document are properties of their registered owners.

Ottobock®
QUALITY FOR LIFE

Malleo Sensa
50S5

Genu Sensa
50S5


Lumbo Sensa
50R10 / 50R11

Genu Sensa
50S5

Manu Sensa
50P13


Sizing Tables

Malleo Sensa 50S5

| | | | | | | | | |
|---|----------------------------------|-------------|-------------------|-------------------|-------------------|-------------------|--------------------|---------------------|
|  | <div>left</div> <div>right</div> | Size | XS | S | M | L | XL | XXL |
| | | cm/ inch | 17-19/ 6.7-7.5 | 19-21/ 7.5-8.3 | 21-23/ 8.3-9.0 | 23-25/ 9.0-9.8 | 25-27/ 9.8-10.6 | 27-29/ 10.6-11.4 |


Colour: grey

Genu Sensa 50K15

| | | | | | | | | | |
|---|------|-------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
|  | f | cm/ inch | 35 - 38/ 13.8-15.0 | 38 - 41/ 15.0-16.1 | 41 - 44/ 16.1-17.3 | 44 - 47/ 17.3-18.5 | 47 - 50/ 18.5-19.7 | 50 - 53/ 19.7-20.9 | 53 - 56/ 20.9-22.0 |
| | Size | | XXS | XS | S | M | L | XL | XXL |
| A= 14 cm 5.5 inch B= 12 cm 4.7 inch | d | cm/ inch | 25 - 28/ 9.8-11.0 | 28 - 31/ 11.0-12.2 | 31 - 34/ 12.2-13.4 | 34 - 37/ 13.4-14.6 | 37 - 40/ 14.6-15.8 | 40 - 43/ 15.8-16.9 | 43 - 46/ 16.9-18.1 |


Colour: grey

Lumbo Sensa 50R10 /
Lumbo Sensa Women 50R11

| | | | | | | | | |
|---|------------------|-------------|---------------------|---------------------|----------------------|-----------------------|-----------------------|-----------------------|
|  | <div>women</div> | Size | XS | S | M | L | XL | XXL |
| | | cm/ inch | 70-80/ 27.6-31.5 | 80-90/ 31.5-35.4 | 90-100/ 35.4-39.4 | 100-110/ 39.4-43.3 | 110-120/ 43.3-47.2 | 120-130/ 47.2-51.2 |


Colour: grey

Epi Sensa 50A7

| | | | | | | | | |
|---|-------------|-------------------|-------------------|-------------------|-------------------|--------------------|---------------------|---------------------|
|  | Size | XXS | XS | S | M | L | XL | XXL |
| | cm/ inch | 17-19/ 6.7-7.5 | 19-21/ 7.5-8.3 | 21-23/ 8.3-9.1 | 23-25/ 9.1-9.8 | 25-27/ 9.8-10.6 | 27-29/ 10.6-11.4 | 29-31/ 11.4-12.2 |

Colour: grey

Manu Sensa 50P13

| | | | | | | | | |
|---|----------------------------------|-------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|
|  | <div>left</div> <div>right</div> | Size | XS | S | M | L | XL | XXL |
| | | cm/ inch | 14-15/ 5.5-5.9 | 15-16/ 5.9-6.3 | 16-17/ 6.3-6.7 | 17-18/ 6.7-7.1 | 18-19/ 7.1-7.5 | 19-20/ 7.5-7.9 |

Colour: grey



- Indications**
- Feeling of instability
 - Irritations (post operative, poist traumatic or chronic)
 - Joint swellings (osteoarthritis or arthritis)
 - Tendomyopathies
 - Ligament insufficiency

- Effects**
- Supports sensorimotor functions of the ankle joint and thus activates muscle stabilization
 - Promotes targeted compression for resorption of edemas, hematomas and effusions
 - Relieves pain due to integrated, anatomical Technogel® pads in the medial and lateral areas
 - Functional improvement of the ankle joint

- Features:**
- Natural movement is not restricted
 - Open mesh section in the heel area for improved wearer comfort
 - Elasticity around edges
 - Can be worn comfortably inside shoe
 - Reduced compression zones (heel & in-step)

Ankle support with pads

- Indications**
- Feeling of instability
 - Irritations (post-operative, post-traumatic, or chronic)
 - Joint swellings (osteoarthritis or arthritis)
 - Knee pain caused by strain (tendomyopathies, insertion ligamentoses)

- Effects**
- Supports sensorimotor functions of the knee joint, relieves pain, and activates muscle stabilization
 - Intermittent compression enables resorption of edemas, hematomas and effusions
 - Relieves and guides the patella due to integrated pad
 - Functional improvement of the knee joint

- Features:**
- Technogel® pad applies pressure to the patella area
 - Lateral spiral springs keep the support from shrinking and provide additional stability
 - Stabilizes the knee
 - The support can be worn on the left or right knee
 - Promotes joint functionality



Patented knit fabric, extremely soft and elastic.

Knee support with pad

- Indications**
- Lower back pain
 - Pseudoradicular lumboischialgia
 - Muscular imbalances in the lumbar spine region
 - Degenerative changes in the lumbar spine region (osteocondrosis, spondylarthrosis, facet syndrome)
 - Irritation of the sacroiliac joint

- Effects**
- Promotes sensorimotor functions and activates muscle stabilization in the lumbar spine region
 - Supports anatomically correct posture
 - Reduces pain and noticeably relieves the lumbar spine
 - The pad has a proprioceptive effect and stabilizes the lumbosacral transition area
 - Technogel® pad improves circulation, stimulates muscle function and reduces muscle tension

- Features:**
- The lumbar support is available in straight style (50R10) and waisted style (50R11)
 - Continuously adjustable closures with pocket grips make the support comfortable to apply



The breathable pad gently massages lumbar spine area, stimulating circulation.

Lumbar support with pad

- Indications**
- Tennis elbow/golfer's elbow
 - Irritations (post-operative, post-traumatic or chronic)
 - Joint swellings (osteoarthritis or arthritis)
 - Tendomyopathies and other soft tissue diseases (chronic bursitis) and injuries

- Effects**
- Provides targeted compression of the elbow region for the resorption of edemas, hematomas and effusions
 - Technogel® pads enhance the micro massage effect on the soft tissue around the joint
 - Relieves pain and improves muscular stabilization

- Features:**
- Technogel® pads in the medial and lateral areas
 - Ring-shaped recesses in the pads serve to relieve the pressure on the epicondyles and contribute to the anatomically correct positioning of the support
 - The support can be worn either on the left or on the right elbow

Elbow support with pads

- Indications**
- Acute and chronic irritations (tendovaginitis, tendomyopathy, minor joint effusions and arthralgias)
 - Post-operative and post-traumatic irritations
 - Sprained wrist
 - Rheumatoid arthritis
 - Carpal tunnel syndrome
 - Degenerative wrist problems (osteoarthritis)
 - Feeling of instability (snapping ulna)

- Effects**
- Compresses the wrist and promotes muscle stabilization
 - Relieves pain and supports functional improvement of the wrist
 - A volar plastic splint ensures optimal positioning of the hand and limits the wrist movement in palmar direction

- Features:**
- Comfort for sensitive wrist bone thanks to open-mesh system
 - The compression reduced section in the knuckle region (styloid process of the ulnar) makes the support comfortable to wear
 - Technogel® pads for optimum compression distribution
 - The volar plastic splint is anatomically pre-shaped and individually adaptable
 - An elastic wrist strap helps to additionally stabilize the wrist
 - Flexible strap can be added for additional wrist joint stabilization

Wrist support with pad