


Rubber Threshold Ramps can be used indoors and outdoors in a variety of locations, whether it's in your home, the office or the garden.

Uses in home may include:

- At front and back entrance doors
- Garage to kitchen
- Kitchen to living room
- Walk-in showers



## Rubber Threshold Ramp

Caution: This product is made from recycled natural rubber which may contain Latex and could cause an allergic reaction in some individuals.

**Pride**<sup>™</sup>  
Lifts & Ramps

**PRIDE**<sup>®</sup>  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST<sup>®</sup> —

www.prideliftsandramps.com • www.pridemobility.com  
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

## Rubber Threshold Ramp



**Pride**<sup>™</sup> Lifts & Ramps

www.prideliftsandramps.com • www.pridemobility.com  
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

# Rubber Threshold Ramp



Rubber Threshold Ramps from Pride™ Lifts and Ramps provide a durable, safe means of navigating uneven household surfaces including small steps and doorways. A hidden dovetail locking system allows for easy installation and incremental ramp height adjustment of up to 3 inches. Modular design makes Pride Rubber Threshold Ramps a practical and affordable home accessibility solution.

- Meets ADA slope requirements
- Constructed of high-impact recycled rubber for durability & long life
- Five total pieces:
  - Increased by 1/2" height increments
  - From 1" to 3" at maximum height
  - 36" of usable space
- Can be used indoors or outdoors
- Slip resistant
- Easy to install
- Maintenance-free
- 850 lbs. weight capacity
- Can be trimmed or notched for a custom fit

## Rubber Threshold Ramps Specifications

RISE	TOTAL LENGTH	USABLE WIDTH	OVERALL WIDTH	TOTAL WEIGHT
1"	10.43"	36"	39"	9.5 lbs.
1.5"	17.33"	36"	41"	21.5 lbs.
2"	24.23"	36"	42.5"	37 lbs.
2.5"	31.13"	36"	44.2"	55 lbs.
3"	38.63"	36"	46"	78 lbs.

The Rubber Threshold Ramp has a maximum incline of 4° (7%), which meets the requirement established by the ADA.



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.