



Technical modifications reserved.

Otto Bock®

QUALITY FOR LIFE

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt/Germany · Phone +49 5527 848-1424 · Fax +49 5527-72330 · export@ottobock.de · www.ottobock.com

© Otto Bock · 646D 439-GB-01-1002

Otto Bock[®]

QUALITY FOR LIFE



Malleo Sprint

Product Information



Malleo Sprint



The number one choice for chronic instability

The new Malleo Sprint ankle orthosis is the number one therapy choice in case of chronic ligament instability in the foot/ankle region. Thanks to the special plastic exoskeleton and strap guide, the Malleo Sprint assures a high degree of stabilisation without having a negative impact on plantar and dorsal flexion. Undesirable eversion and inversion of the foot are limited.

The orthosis structure, which adapts perfectly to the anatomical shape of the body, assures targeted, stable support and the limitation of movement when it is closed. Breathable materials and excellent padding enhance the level of wearer comfort.

Malleo Sprint

Features and Benefits

High level of wearer comfort thanks to an outstanding anatomical fit

An anatomical protrusion in the malleoli area assures enhanced wearer comfort. The reduced thickness of the base material under the foot combined with highly comfortable padding prevent unwanted pressure points.



High degree of stability

The combination of a plastic skeleton and strap guide assures a high degree of stability and prevents unwanted movements.



Optimum hook closure system

The integrated hook guide assures optimum support for the laces and allows the orthosis to be put on and taken off quickly.




Easy entry

A grip loop helps in putting the orthosis on quickly and comfortably.



Indications

- Chronic ankle joint instability
- Functional treatment of ligament ruptures
- Severe ankle distortions



Order no.	Size	Circumference (cm)
50S3=XS-1	XS-1	approx. 25.5 – 28.0
50S3-1=S	S-1	approx. 28.0 – 30.5
50S3=M-1	M-1	approx. 30.5 – 33.0
50S3=L-1	L-1	approx. 33.0 – 35.5
50S3=XL-1	XL-1	approx. 35.5 – 38.0

The orthosis is universal for the left or right side.