

Otto Bock®

QUALITY FOR LIFE



Manu Immobil Manu Immobil Long

For occupational or
post-operative Treatment

The 50P10 Manu Immobil and 50P11 Manu Immobil Long wrist positioning orthoses has been designed for immobilization of the hand and wrist in a functional position. The 50P11 Manu Immobil Long offers additional fixation of the finger joints. Both orthoses can be used for post-operative and post-traumatic immobilization as well as for occupational therapy.

Use of the special three-dimensional SpaceTex® spacer fabric ensures high comfort for the patient. SpaceTex® is breathable and has excellent padding qualities.

The orthoses are lightweight and – since the pad and straps are washable – offer good hygienic properties.

Otto Bock HealthCare GmbH

Max-Näder-Straße 15 · 37115 Duderstadt · Phone +49 5527 848-0
healthcare@ottobock.de · www.ottobock.com



Manu Immobil 50P10

Manu Immobil Long 50P11

Indications

- After fractures and luxations of the hand, the wrist (applies to 50P10 and 50P11) and of the finger joints (only applies to 50P11)
- After wrist, carpal and metacarpal surgery
- For treatment of acute irritations (such as carpal tunnel syndrome), osteoarthritis of the wrist or severe distortions (sprains)

The indication must be determined by the physician.

Effects

The hand or wrist and finger joints (only applies to 50P11) are immobilized in a functional position in order to prevent harmful movements.

Features

- Good immobilization in a functional position
- The integrated steel bar offers an individual adjustment in the area of the wrist and metacarpophalangeal joints
- Excellent padding qualities thanks to the three-dimensional SpaceTex® spacer fabric:
 - Air permeability
 - Good adjustability
 - High comfort
- Washable pads and straps
- Lightweight




50P10

50P11

Fax Order Form: +49 5527 72330

Sender/Receiver of goods

| Art. no. | Size | Hand length A | Quantity |
|----------|-----------|----------------------------|----------|
| 50P10 | universal | universal | Pc. |
| 50P11=M | M | up to 18 cm (7.09 inch) | Pc. |
| 50P11=L | L | from 18 cm (7.09 inch) | Pc. |



The orthoses can be used for left or right.

Customer ref.

Customer no.