

Pride Lift Chairs

Experience the Luxury

The LC-108 is now part of the Specialty Collection of Lift Chairs from Pride. Here are all of the Lift Chairs in Pride's Specialty Collection.

SPECIALTY COLLECTION

SPECIALTY COLLECTION MODELS

SPECIAL FEATURES



SPECIAL VALUE



SPECIFICATIONS

	SPECIAL FEATURES										SPECIAL VALUE			
	GL-310	LC-525PW	LC-525S	LC-525M	LC-525L	LL-570T	LL-571	LL-660	LL-670	LL-805	LC-102 ³	LC-106	LC-107	LC-108
Seat-To-Floor:	20"	19"	19"	20"	21"	20"	21"	19"	20"	19"	16"	18"	19"	19"
Seat Depth:	20"	19"	18"	20"	22"	20"	22"	20"	21"	20"	17"	19"	19"	19"
Seat Width:	21"	24"	18"	20"	22"	20"	22"	18"	20"	20"	16"	20"	20"	20"
Top of Back to Seat:	24.5"	25.5"	25.5"	26.5"	27.5"	29.5"	27"	25"	27"	27"	23.5"	27.5"	27.5"	29"
Distance From Wall:	18"	17"	17"	17"	17.5"	20.5"	16"	25"	27"	4"	16"	18"	26"	24"
Overall Chair Width:	30"	39"	35"	35"	37"	35.5"	38"	34.5"	38"	35"	25.5"	32.75"	35"	31.5"
Chair Weight:	115 lbs.	131 lbs.	122 lbs.	123 lbs.	127 lbs.	128 lbs.	133 lbs.	127 lbs.	134 lbs.	118 lbs.	106 lbs.	96 lbs.	110 lbs.	136 lbs.
Weight Capacity:	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	300 lbs.	250 lbs.	325 lbs.	325 lbs.	375 lbs.
Hand Control:	Easy touch	6-button infinite position	6-button infinite position	6-button infinite position	6-button infinite position	Easy touch	Easy touch	Infinite-position toggle	Infinite-position toggle	Easy touch	Toggle	Easy touch	6-button infinite position	6-button infinite position
Warranty:	7-year prorated	2-year limited	2-year limited	2-year limited	2-year limited	7-year prorated	7-year prorated	7-year prorated	7-year prorated	4-year limited	2-year limited	1-year limited	2-year limited	2-year limited

FABRIC SELECTION



BENEFITS

- Engineered furniture grade laminate/hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head and arm covers (head cover only on GL-310)
- Fire retardant (California specifications; bulletin 117)
- Standard, large dual pockets
- Removable back - no tools required (LC-108 requires tools)
- Integrated battery backup in case of power failure²
- Transformer features self-diagnostic electronics with 12' cord
- New "zero-gravity" position with active seat box for maximum comfort (only available on LC-525 Series)

OPTIONS*

- Heat & Massage
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Bed option with mattress, fitted sheet & storage bag (LL-670 only)
- Leather/vinyl wipes
- Protective pad
- Non-skid leg levelers

* Not all options are available for all models. Check order form for details. Additional charges apply for all options.



Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

**5-year warranty.



www.pridemobility.com

St. Catharines, Ontario • 1-888-570-1113

www.pridemobility.com

St. Catharines, Ontario • 1-888-570-1113

Pride
Mobility Products Co.

Engineered for comfort and durability, the LC-108 Infinite-Position Lift Chair from Pride® reclines to a full Trendelenburg position that elevates the legs above the head for maximum comfort and relaxation. The stylish LC-108 is available in a range of eye-catching fabric selections sure to enhance your home's décor.

NEW! **LC-108**

Pride Lift Chairs | Experience the Luxury

Designed for ultimate comfort and adjustability to meet individual needs

The LC-108 is equipped with adjustable and removable components to allow for a custom fit and easier care and maintenance of the product.

Adjustable, removable back cushion

- **The individual pillows of the back have zippers** to access the fiber fill stuffing. The stuffing can be added or taken away to customize the fit and comfort of the back cushions. This is especially helpful with accommodating a user that has a kyphosis. For users with longer legs, stuffing can be removed from the back cushions to increase seat depth, offering better overall comfort and support for the legs.
- **The individually tethered pillows** of the back allow them to move with the client when going into or coming out of recline to reduce shearing and shifting.
- **The standard wing back design** makes it possible to insert additional support along the side to address issues with hemiplegia or stroke. The wings also work as a safeguard to keep users from falling beyond the sides of the chair.



WING BACK DESIGN



REMOVABLE / ADJUSTABLE
BACK CUSHION



ZIPPERS TO ACCESS
FIBER FILL STUFFING



INDIVIDUALLY TETHERED
PILLOWS

Adjustable, removable chaise pad

- **The chaise pad is fully removable** and the foam fill can be removed. This will allow for hand washing and hang drying of the seat cushion should it become soiled. The foam fill must be removed during washing. Moisture can cause damage to the foam.
- **Access to the seat cushion** also gives the ability to insert seating options to address pressure or comfort concerns.
- **The full chaise pad** provides a continuous support surface for the entire length of the seating surface, reducing pressure points and allowing for maximum comfort and relaxation.
- **Seating wedges or components** can be put between the seat box and chaise pad to correct postural seating concerns.



FULLY RECLINED WITH
BACK AND CHAISE CUSHIONS



REMOVABLE CHAISE PAD



BACK AND CHAISE
CUSHIONS REMOVED

Additional Features

- **Both the seat and back cushions are replaceable** as parts should it be necessary due to wear and tear, increasing the life of the chair. (Note: Over time fabric fades and with different dye lots of fabric, replacement parts can sometimes differ slightly in color).
- **Push-button hand control standard** that allows for independent positioning of the leg rest and back angle. This control is also equipped with a button that, when depressed, will lower the legs, raise the back and then lift the chair to the transfer position all in one continuous motion. The hand control has a quick disconnect feature in the side pocket to disable the functions of the control.
- **Low-shear recline back standard.** This means that there is minimal shifting when going into recline and less forward force on the client when coming out of recline. This reduces shearing and maintains better positioning and comfort for the user.
- **Fully landscaped arms** provide a continuous, even arm support surface with no ridges or transitions in height. This offers better comfort and arm support.
- **The standard footrest extension** gives more support to the feet, reducing pressure on the back of the calves and increasing relaxation to the feet. This often helps address cold feet issues.
- **Trendelenburg or zero gravity positions** are achieved when the chair is fully reclined and the feet are fully elevated. In this position, the seat button can be used to shift the entire seating system to raise the head and lower the legs while maintaining the seat to back angle. This allows for greater range of seat positions to personalize comfort and select the best napping position.
- **Standard arm and head rest covers** are supplied with the chair. This allows for easier care and maintenance of the chair and increases fabric protection.
- **Optional head rest and lumbar roll** can be added to the chair to further enhance comfort and positioning.



TRENDELENBURG / ZERO GRAVITY
POSITION



SHOWN WITH STANDARD
HEAD AND ARM COVERS



PUSH BUTTON HAND CONTROL